



Transportation Bag STD for Helicopter Rescue Bag Set



For the operation-ready transport of a helicopter rescue bag, a vacuum mattress, a vacuum pump and additional accessories (e.g. RBS).

The transport bag is equipped with a fully stowable **shoulder strap system** and **straps** for hanging on the load hook.

This means that the complete system can also be carried over longer distances if necessary.

The rescue bag is folded twice and placed in the transport bag with the vacuum mattress inserted. The bag can be opened quickly and clearly thanks to the **all-round, robust zipper**. The rescue bag is removed, unfolded and is immediately ready for use.

Made from **tear-resistant, water- and wind resistant special Perlon fabric**. With all-round zip and storage pocket for the belt system. Two detachable **compression straps** allow the pack size to be further reduced. The back area is additionally reinforced with Cordura.

Sufficient space for:

- Helicopter Rescue Bag
- Vacuum Mattress
- Vacuum Pump
- RBS
- Face Protection Shield
- Small Accessories



Shoulder-Strap System

in use

stored-away



P/N	93227	
Dimensions	Width:	600 mm
	Height:	650 mm
	Depth:	300 mm
Weight	1,40 kg	
Comments	Also available as a SMALL transportation bag to stow only the rescue bag (# 93395).	