

# TECH - SPEC

## MODEL 350 EXTREME RESCUE TOBOGGAN

### Standard Features

- ◆ Four fixed position removable handles
- ◆ Sure grip finish with "T" handles
- ◆ Chain Brake Release
- ◆ Three patient restraint straps
- ◆ Weighs 55 lbs

Length 93", Width 26"

### Options

- ◆ Rated Belay Points (*Independently tested*)
- ◆ FastPax Restraint System
- ◆ QuickTab 1.5" patient restraint straps
- ◆ Carbon fiber handles (*8lbs - set of four*)
- ◆ Two piece version
- ◆ Two piece Carbon fiber version (*33.5lbs*)



Copyright 2017 All Rights Reserved - Cascade Rescue Company  
1808 Industrial Drive, Sandpoint, ID 83864  
Phone: (208) 263-2484 Fax: (208) 263-255-7460 [www.cascade-rescue.com](http://www.cascade-rescue.com)