




## CARRYING STRAPS TYRAL

Cross carrying straps for a more **efficient distribution of the weight** to be transported when removing an injured person using the mountain stretcher TYRAL.

The **padded shoulder pads** and the adjustable and also in the dark well visible holding loops enable a safe and comfortable carrying of the GT TYRAL. The straps can be quickly and easily adjusted by means of the robust buckles, depending on body size and trail conditions.

They can **simply be hooked into the notch** between the tube and the oval handle. This makes it easy to change the rescuer carrying the stretcher at any time.

Additionally available: **loop with cobra buckle** for quick fixation (#93259).



Art. No.	<b>93256</b>
Dimensions	
Weight	0,3 kg
Comments	Supplied as a set.