



Tragekraxe TYRmove Akja version



Struts and belts for für optimale fixation of half- Akja.



Two pair of poles can be fixed within the straps. Therefore no shifting of the load is possible.



Comfortable, energy-saving carrying case for high alpine use.

Easily attach and carry rescue stretcher, divisible akja, vacuum mattress and the like.

Special connection with strut and fixing belt. The pair of drawbars can be attached to the fixation strap using special strap eyelets.

Your own backpack with the man's equipment can also be fixed.

The shoulder strap and hip belt of the Kraxe are well padded and serve as protection against pressure points. The V2 carrying system is dimensionally stable and durable.

No shifting of the load while walking - especially important in inaccessible terrain. The load is fastened using three fixing straps with a stable clamp. Very robust and easy to use system.

Sturdy aluminum frame and comfortable carrying system enable versatile use - load capacity up to 50 kg.

A 25 cm deep arm with stand serves as the loading area. The stand enables a secure stand - even on uneven terrain.

This also offers additional mounting options.

Scope of delivery: two pieces of carrying case TYRmove incl. two struts and fixing straps per piece.

Art. Nr.	93800	
Dimensions	Length:	780mm
	Width:	350mm
	Stage space:	250mm
Weight	2x4kg	
Comments	For each half of an Akja divisible two pieces of carrying bag are necessary	