



## TRANSPORTBAG TYRAL CAR

For **protected transport** and **safe storage** of the TYRAL mountain stretcher.

The transport bag TYRAL is made of **waterproof tarpaulin fabric**, which keeps the mountain stretcher dry and clean even during transport on the car roof.

There are two zippers on the opening flap. These can be opened over the entire length.

**Easy and quick loading / unloading** of the stretcher is thus enabled.

Circumferentially there are belt loops. These are provided for a safe fastening.

For fastening and transport on public roads, a lashing device must be used, which complies with the **ÖNORM EN 12195-2** standard.

When fastening, it must be ensured that not only the transport bag alone, but also the stretcher in the interior is tensioned down. The transport bag serves exclusively to protect the stretcher during transport.



Art. No.	<b>93740</b>
Dimensions	Length: 2020 mm Width: 650 / 300 mm
Weight	3,5 kg
Comments	Without fastening straps.