



INJURED PROTECTION BAG LITE



Light, **weight-optimized** protection bag for injuries in a compact design.

Suitable for **quick, short rescues** from easy terrain.

Hood in the head area for additional **weather protection**.

The **foot section** with a separate opening allows flexible adaptation to voluminous ski boots or rails.

Supplied with the **TYRbag 4L**, making it easier to **transport** the protective bag.



Art. Nr.	93086
Dimension	Length: 2050 mm Width: 450 mm
Weight	1,5 kg
Comment	