



Rescue Trailer TYROMONT for snow-mobiles and quads - 4X4wheels



High-grade, completely equipped rescue trailer for snow-mobiles and quads.

Optimum protection of the injured against rain, snow-spray, wind and cold.

Maximum safety by two rollover bars and fixation belts with aluminum buckles.

Integrated, extractable stretcher with adjustable reclining area with thermo-insulation. Stretcher equipped with diagonal fixation belts for the upper part of the body and transverse belts for the hip and legs.

Additionally equipped with two **ground anchors** for fixation of the stretcher on snow or frosty underground.

Guiding profiles and **quick-locks** for easy loading of the injured on the stretcher.

Additional **guiding profiles** for optional loading of an Akja or Canadier rescue sledge (please specify your demand during ordering).

Cover manufactured from water- and wind-resistant PES fibre with large, cold-resistant clear-areas.

Cover could be opened by a **large, framed zipper** to easily handle the injured person.

Stable running properties through 4 Balloon wheels.

Dampened suspension for a smooth ride of the injured.

| | |
|--------------------|---|
| Art. No. | 93158 |
| Dimension | Total : 2.250 x 790mm Length trailer: 2.850 mm Width incl. wheels: 1.500 mm Height incl. cover: 1.470 mm |
| Weight | Trailer incl. rescue set: 105 kg |
| Also available as: | <ul style="list-style-type: none"> • Use with Canadier and injured stretcher • After-market equipment without trailer |