



Rescue Trailer TYROMONT for snow-mobiles and quads



High-grade, completely equipped rescue trailer for snow-mobiles and quads.

Optimum protection of the injured against snow-spray, wind and cold.

Maximum safety by two rollover bars and fixation belts with aluminum buckles.

Integrated, extractable stretcher with adjustable reclining area with thermo-insulation. Stretcher equipped with diagonal fixation belts for the upper part of the body and transverse belts for the hip and legs.

Additionally equipped with two **ground anchors** for fixation of the stretcher on snow or frosty underground.



Injured stretcher fixed on trailer.



Guiding runners with fins.



Ground anchors for fixation of stretcher on a slope.

Guiding profiles and **quick-locks** for easy loading of the injured on the stretcher.

Additional **guiding profiles** for optional loading of an Akja or Canadier rescue sledge (please specify your demand during ordering).

Cover manufactured from water- and wind-resistant PES fibre with large, cold-resistant clear-areas.

Cover could be opened by a **large, framed zipper** to easily handle the injured person.

Stable running properties through wide track, guiding runners and additional four fins. **Dampened suspension** for a smooth ride of the injured.

Art. No.	93152
Dimension	Total : 2.250 x 790mm
	Length trailer: 2.850 mm
	Width incl. runners: 1.080 mm
	Height incl. cover: 1.450 mm
Weight	Trailer incl. rescue set: 105 kg
Also available as:	<ul style="list-style-type: none"> Use with Canadier and injured stretcher After-market equipment without trailer