



COMBI ARM - LEGSPLINT

Multi-use combination splint for immobilization. This easy-to-combine **arm/leg splint** has a modular design. Multiple injury patterns can be immobilized with this splint.

The armsplint is attached to the legsplint with Velcro and thus serves as a large leg splint. Due to the two halves of the splint, more precise shaping is safely possible even with only two rescuers.

Arm section can also be used as a leg splint for children, but also for forearm injuries in adults.

Leg section is large enough to adequately immobilize ankle or knee injuries, for example.

Automatic, self-shutting valve does not have to be closed after vacuumizing the mattress. Even if the pump hose slips off, there is no loss of vacuum. **Two separate valves for suction.**



Art. No.	81354	
Dimensions	Length:	1100 mm
	Width:	660 mm
Weight	1,35 kg	
Comments		