

## TYROMONT Appine Rescue Equipment

## Akja 2000 DIVISIBLE







The classic Akja - Rescue Sled in **divisible version** in length of 200 cm.

Designed for **2-rescuer** operation in **all kind of snow** - covered terrain, groomed or ungroomed. **Divisible in two halves** for transport on the back in unexploited areas.

**Long lasting durability:** Manufactured completely from high-strength, corrosion-resistant aluminum alloy for many years of reliable, maintenance-free action.

Solid and **quick dividing - mechanism** provides a safe hold while riding and quick parting in the middle if necessary for transport.

**Stable running characteristics**: Reinforced steel runners with especially hardened fins provide stable running even in very icy conditions.

**Two pairs of poles** (for two rescuers), can be detached easily by bayonett-lockings. The grips of the handles are high enough for tall rescuers. With integrated stabilization profiles for increased bending strength.

Optionally, divisible handles are available.

**Bayonett-lockings** are wear-free, 100% icing-safe and easily to handle.

**Carrying grips** on both ends of the sledge for comfortable lifting and carrying of the sled.

**Four fixation belts with aluminum buckles** for safe fixation of the injured; resistant against breaking and icing, can be operated with skiing-gloves.

Art. No.	93032	
Dimension	Length: Width: Length divided:	2.000 mm 560 mm 1.000 mm
Weight	27,0 kg (incl. poles)	
Comments	Divisible handles optionally available.	