



## 4f Adjustable Suspension TYRAL



**Adjustable suspension** for TYRAL mountain stretcher.

Adjustable suspension so that the patient lies as horizontally as possible in the stretcher.

Loops in the head and leg area can be adjusted from 30cm to 150cm.

The four adjustable slings allow the stretcher to be the stretcher can always be ideally aligned depending on the terrain and weight shift.

The suspension is attached to the 4 attachment points at the head/foot area using twistlock carabiners.

The rope is attached by means of a rigger plate and carabiner.

All parts certified according to EN 354, EN 566 and EN 795 -B.



|            |                             |           |
|------------|-----------------------------|-----------|
| Art. No.   | <b>93512</b>                |           |
| Dimensions | Length head:                | 30-150 mm |
|            | Length foot:                | 30-150mm  |
| Weight     | 1,6 kg                      |           |
| Comments   | Supplied with TYRbag pouch. |           |