



Salvage transport bag 7- Terra



The **salvage transport bag 7-Terra** has been designed for the transport of fatally injured persons in the alpine terrain.

It is characterized by its light weight and easy handling. The **12 carrying loops** make it possible to carry the casualty safely in pathless terrain.

The large zipper on three sides allows easy opening and closing of the transport bag at any time. On the front side there is an **A5 document pocket** with a large viewing window for storing important information.

The **welded seams** inside the bag prevent water and other liquids from escaping unintentionally.

The use of tarpaulin fabric allows **easy and quick cleaning** with water and a light soap suds without long drying time.



Art. No.	93090	
Dimensions:	Length:	2020 mm
	Width:	500 mm
Weight:	3,0 kg	
Comments:	Supplied with transport bag.	