



## Thermo Rescue Bag - TRB



Newly developed rescue bag for the ideal transport of injured in alpine terrain.

**High-quality rescue bag** for use with Akja, Canadier, Mountain Stretcher or any other type of rescue device.

Ergonomic design, **high-strength, water-resistant nylon-fabric**, with strong wide area Velcro fastener for adjustable closing of the bag.

**Two additional zippers** at the upper part of the bag allow a fast access to the injured person (e.g. for an infusion).

**Foot-compartment with two zippers** can be opened separately and provides enough space for large ski boots and/or splints.

**Thermo insulation layers** at both sides and the bottom that can be removed for cleaning purposes.

**Six tight carrying loops** for comfortable carrying of the injured.

Very **compact foldable**, little weight of 2,0 kg.

Without additional dampening mat at the bottom; therefore, optimum application in **combination with a vacuum mattress**.

Suitable for all TYROMONT rescue sleds, mountain stretchers and other models.



Art. No.	<b>93070</b>	
Dimension	Length:	2.050 mm
	Width:	550 mm
Weight	2,0 kg	
Comments	Ideal combination with a vacuum mattress.	