



Mountain Stretcher STEELIGHT

Extremely rugged and all-purpose mountain stretcher for tough, alpine rescuing operations.



Solid mountain stretcher for the transport of injured persons in all kind of terrain, for roping up or down, dragging over rock or ice, carrying and riding on the single wheel - wherever a good protection of the injured and a rugged device are required.

New ideas and well-proven design

The proved, circle-sector shaped frame design is the best choice for easy handling in all kind of terrain. The flat lying area for the injured is totally protected from side-impacts by a rail that surrounds the whole stretcher. This rail also allows to hold the stretcher at any position. The two base runners are shaped in a circle-sector for optimum handling and retracted to the center to provide more space for the legs of the rescuer when roping up or down in steep rock. For transport reasons the mountain stretcher can be quickly divided in the middle with Camlock quick locks.

Adjustable handles for more convenience

The four handles can be quickly and individually adjusted in 12 different positions to adapt to different carrying situations. The adjustment of the handles are done by arresting bolts. The handles can be fully removed.

Protection and best care for the injured

The absolutely horizontal bedding area for the injured is made from extra-light but extremely impact-resistant aluminium sandwich material, covered with shock-absorbing foam. It is 2.000 mm long, the width is 500 mm in the hip, 530 mm in the shoulder and 380 mm in the head and foot area, providing enough space even for tall patients. It allows unrestricted use of a vacuum mattress.

Four fixation belts with aluminium buckles hold the injured safe in every situation.

Extended field of application. The mountain stretcher light comes complete with wheel struts and a high-strength 20 inch aluminium wheel with a tubeless tire. The new developed wheel strut EVO could be easily attached on the underframe of the stretcher and is fastened with two ball-lock bolts. Foldable to reduce size during transport.



Art. No.	93240	
Dimensions	Length:	2.000 mm
	Width:	530 mm
Weight	16,0kg	(stretcher w/o handles)
	5,0 kg	(handles)
	3,0 kg	(wheel with struts)