



Vacuum Mattress

ULTRALIGHT

with Fixation Belts and Carrying Handles

This new type of vacuum mattress offers several **advantages** over the previous vacuum mattresses:

Light weight only **4.7 kg**.

New multi-chamber system - mattress maintains its shape. No more stroking necessary.

New **removable shoulder straps** and carrying straps.

Can be **rolled** up and easily transported with a shoulder strap.

The new **multi-chamber system (120 chambers)** - the layers of the outer cover are welded together point by point, creating chambers between the points that hold the inner beads in their correct position and are evenly distributed.

This enables **quick, safe and easy use** in all situations.

Automatic, self-shutting valve does not have to be closed after vacuumizing the mattress. Even if the pump hose slips off, there is no loss of vacuum. The valve can be operated with any commercially available vacuum pump using an adapter.

Temperature resistance: -40°C to +70°C



Art. No.	81370	
Dimensions	Length:	1.950 mm
	Width:	900 mm
Weight	4,7 kg	
Comments	Valve at shoulder section.	