

# TECH - SPEC

## CRUISER SS RESCUE TOBOGGAN



Length 93", Width 22"



### Standard Features

- ◆ Variable Height Handles
- ◆ Integrated Spray Guard
- ◆ On the fly Chain Brake Release
- ◆ Sure Grip knurled handles with "T" handles
- ◆ Insulated, padded interior
- ◆ Snow Slick bottom for increased performance
- ◆ Advanced Fin design
- ◆ Stainless steel hardware throughout
- ◆ Weighs 52 lbs

### Options

- ◆ Rated Belay Points (*Independantly tested*)
- ◆ FastPax Restraint System
- ◆ QuickTab 1.5" patient restraint straps
- ◆ Slide Arrest Anchors
- ◆ Stainless steel Bumper Handles
- ◆ Advanced Fin design
- ◆ Stainless steel hardware throughout

