

TECH - SPEC

CASCADE TOBOGGAN MODEL 100 "THE LEGEND" RESCUE TOBOGGAN



Standard Features

- ◆ 30 or 50 Inch Fins
- ◆ Multi-Layer Composite Shell
- ◆ Stainless Steel Perimeter Frame and Handle Mounting Components
- ◆ Stainless Steel Handles
- ◆ 6061 Aircraft Aluminum Handle Locks
- ◆ Three Patient Restraint Straps
- ◆ Quick Detach Tail Rope
- ◆ Chain Brake
- ◆ Four Carry Handles
- ◆ Crossover Rope
- ◆ Non-Slip Interior

Options

- ◆ Bow Guards - Reduce wear and tear and increase the life of your equipment
- ◆ Sled Loader - No need to "lap carry". Safe, stable and easy to use.
- ◆ Snowmobile Tow Bar - Standard or Stow -N- Go
- ◆ Chain Brake Release, Gripper Chain, Belay Points and Snow Slick
- ◆ Patient Protector, Pad Stretcher,

Specifications

Model	The Legend, Model 100 w/30" Fins	The Legend, Model 100 w/50" Fins
Item	CTM100-30	CTM100-50
Length-Handles Folded	91 in. (231cm)	91 in. (231cm)
Length-Handles Ext.	150 in. (381cm)	150 in. (381cm)
Width	21.5 in. OD, 18.5 ID (54/47cm)	21.5 in. OD, 18.5 ID (54/47cm)
Depth	5 in. (12.7 cm)	6 in. (12.7 cm)
Weight	52 (23.6kg)	52 (23.6kg)

