



## Vacuum Mattress 60 with Cover and Carrying Handles



**New vacuum mattress:** designed for the special demands of mountain and air rescuing.

**Compact shape:** provides extremely high side and length stability at very compact measures. Dimensions exactly adjusted to TYROMONT rescue sledges, mountain stretchers and helicopter rescue bags.

**Revolutionary 40-chamber system:** Holds the filling material always at its place. That's why the vacuum mattress can be instantly applied in most difficult terrain.

**High-grade filling material** and the 40-chamber system allow an extremely low profile of the vacuum mattress at maximum stability in vacuumized state. Consequently, the mattress can be shaped exactly to the body of the injured and thus provides optimum immobilization.

**Stable inside-hull** made from exceptionally resistant, fibre reinforced PU-material. Durable welds thanks to CAM welding-procedure.

**Low packing volume** because of the little height of 50 mm of the mattress.

**Rugged outside-cover** provides extra protection for the mattress and is equipped with **eight carrying belts**. For cleaning and disinfection, the cover could easily be removed.

**Automatic self-shutting valve** does not need to be closed after vacuumizing the mattress. Also, the vacuum remains unbiased if the hose of the pump should slip of.

**Temperature resistance:** - 40°C to +70°C

**Manufactured according ISO 9001 and CE.**

Art. No.	<b>81354</b>	
Dimensions	Length:	1.950 mm
	Width:	600 mm
Weight	4,9 kg	
Comments	Valve at foot section.	