



SUSPENSION PRO TYRAL




Suspension Pro TYRAL mountain stretcher.

Sling leg area 108 cm; Orange 21mm
Loop head area 100 cm; Yellow 21mm

Attachment of the stretcher using a Twistlock carabiner
Attachment to the rope on 2 rigger plates and carabiners.

The suspension can be converted from a horizontal to a vertical suspension without tools.
Simply unthread the two yellow ribbon loops from the rigging plate and use them for vertical lowering.



Art. Nr.	93257	
Dimension	Length Head:	100 mm
	Length Feetarea:	108 mm
Weight	total:	0,8 kg
Comment	Delivered with TYRbag	